

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
<i>ABL1</i>	✓	-	✓
<i>ABL2</i>	✓	-	-
<i>ACVR1</i>	✓	-	-
<i>ACVR1B</i>	✓	-	-
<i>AKT1</i>	✓	-	-
<i>AKT2</i>	✓	✓	-
<i>AKT3</i>	✓	-	✓
<i>ALK</i>	✓	✓	✓
<i>ALOX12B</i>	✓	-	-
<i>ANKRD11</i>	✓	-	-
<i>ANKRD26</i>	✓	-	-
<i>APC</i>	✓	-	-
<i>AR</i>	✓	✓	✓
<i>ARAF</i>	✓	-	-
<i>ARFRP1</i>	✓	-	-
<i>ARID1A</i>	✓	-	-
<i>ARID1B</i>	✓	-	-
<i>ARID2</i>	✓	-	-
<i>ARID5B</i>	✓	-	-
<i>ASXL1</i>	✓	-	-
<i>ASXL2</i>	✓	-	-
<i>ATM</i>	✓	✓	-
<i>ATR</i>	✓	-	-
<i>ATRX</i>	✓	-	-
<i>AURKA</i>	✓	-	-
<i>AURKB</i>	✓	-	-
<i>AXIN1</i>	✓	-	-
<i>AXIN2</i>	✓	-	-
<i>AXL</i>	✓	-	✓
<i>B2M</i>	✓	-	-
<i>BAP1</i>	✓	-	-
<i>BARD1</i>	✓	-	-
<i>BBC3</i>	✓	-	-
<i>BCL10</i>	✓	-	-
<i>BCL2</i>	✓	-	✓
<i>BCL2L1</i>	✓	-	-
<i>BCL2L11</i>	✓	-	-
<i>BCL2L2</i>	✓	-	-
<i>BCL6</i>	✓	-	-
<i>BCOR</i>	✓	-	-
<i>BCORL1</i>	✓	-	-
<i>BCR</i>	✓	-	-
<i>BIRC3</i>	✓	-	-

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
<i>BLM</i>	✓	-	-
<i>BMPR1A</i>	✓	-	-
<i>BRAF</i>	✓	✓	✓
<i>BRCA1</i>	✓	✓	✓
<i>BRCA2</i>	✓	✓	✓
<i>BRD4</i>	✓	-	-
<i>BRIP1</i>	✓	-	-
<i>BTG1</i>	✓	-	-
<i>BTK</i>	✓	-	-
<i>C11orf30</i>	✓	-	-
<i>CALR</i>	✓	-	-
<i>CARD11</i>	✓	-	-
<i>CASP8</i>	✓	-	-
<i>CBFB</i>	✓	-	-
<i>CBL</i>	✓	-	-
<i>CCND1</i>	✓	✓	-
<i>CCND2</i>	✓	-	-
<i>CCND3</i>	✓	✓	-
<i>CCNE1</i>	✓	✓	-
<i>CD274</i>	✓	-	-
<i>CD276</i>	✓	-	-
<i>CD74</i>	✓	-	-
<i>CD79A</i>	✓	-	-
<i>CD79B</i>	✓	-	-
<i>CDC73</i>	✓	-	-
<i>CDH1</i>	✓	-	-
<i>CDK12</i>	✓	-	-
<i>CDK4</i>	✓	✓	✓
<i>CDK6</i>	✓	✓	-
<i>CDK8</i>	✓	-	-
<i>CDKN1A</i>	✓	-	-
<i>CDKN1B</i>	✓	-	-
<i>CDKN2A</i>	✓	-	-
<i>CDKN2B</i>	✓	-	-
<i>CDKN2C</i>	✓	-	-
<i>CEBPA</i>	✓	-	-
<i>CENPA</i>	✓	-	-
<i>CHD2</i>	✓	-	-
<i>CHD4</i>	✓	-	-
<i>CHEK1</i>	✓	✓	-
<i>CHEK2</i>	✓	✓	-
<i>CIC</i>	✓	-	-
<i>CREBBP</i>	✓	-	-

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
CRKL	✓	-	-
CRLF2	✓	-	-
CSF1R	✓	-	✓
CSF3R	✓	-	-
CSNK1A1	✓	-	-
CTCF	✓	-	-
CTLA4	✓	-	-
CTNNA1	✓	-	-
CTNNB1	✓	-	-
CUL3	✓	-	-
CUX1	✓	-	-
CXCR4	✓	-	-
CYLD	✓	-	-
DAXX	✓	-	-
DCUN1D1	✓	-	-
DDR2	✓	-	-
DDX41	✓	-	-
DHX15	✓	-	-
DICER1	✓	-	-
DIS3	✓	-	-
DNAJB1	✓	-	-
DNMT1	✓	-	-
DNMT3A	✓	-	-
DNMT3B	✓	-	-
DOT1L	✓	-	-
E2F3	✓	-	-
EED	✓	-	-
EGFL7	✓	-	-
EGFR	✓	✓	✓
EIF1AX	✓	-	-
EIF4A2	✓	-	-
EIF4E	✓	-	-
EML4	✓	-	✓
EP300	✓	-	-
EPCAM	✓	-	-
EPHA3	✓	-	-
EPHA5	✓	-	-
EPHA7	✓	-	-
EPHB1	✓	-	-
ERBB2	✓	✓	✓
ERBB3	✓	✓	-
ERBB4	✓	-	-
ERCC1	✓	✓	-

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
<i>ERCC2</i>	✓	✓	-
<i>ERCC3</i>	✓	-	-
<i>ERCC4</i>	✓	-	-
<i>ERCC5</i>	✓	-	-
<i>ERG</i>	✓	-	✓
<i>ERRF1</i>	✓	-	-
<i>ESR1</i>	✓	✓	✓
<i>ETS1</i>	✓	-	✓
<i>ETV1</i>	✓	-	✓
<i>ETV4</i>	✓	-	✓
<i>ETV5</i>	✓	-	✓
<i>ETV6</i>	✓	-	-
<i>EWSR1</i>	✓	-	✓
<i>EZH2</i>	✓	-	-
<i>FAM123B</i>	✓	-	-
<i>FAM175A</i>	✓	-	-
<i>FAM46C</i>	✓	-	-
<i>FANCA</i>	✓	-	-
<i>FANCC</i>	✓	-	-
<i>FANCD2</i>	✓	-	-
<i>FANCE</i>	✓	-	-
<i>FANCF</i>	✓	-	-
<i>FANCG</i>	✓	-	-
<i>FANCI</i>	✓	-	-
<i>FANCL</i>	✓	-	-
<i>FAS</i>	✓	-	-
<i>FAT1</i>	✓	-	-
<i>FBXW7</i>	✓	-	-
<i>FGF1</i>	✓	✓	-
<i>FGF10</i>	✓	✓	-
<i>FGF14</i>	✓	✓	-
<i>FGF19</i>	✓	✓	-
<i>FGF2</i>	✓	✓	-
<i>FGF23</i>	✓	✓	-
<i>FGF3</i>	✓	✓	-
<i>FGF4</i>	✓	✓	-
<i>FGF5</i>	✓	✓	-
<i>FGF6</i>	✓	✓	-
<i>FGF7</i>	✓	✓	-
<i>FGF8</i>	✓	✓	-
<i>FGF9</i>	✓	✓	-
<i>FGFR1</i>	✓	✓	✓
<i>FGFR2</i>	✓	✓	✓

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
<i>FGFR3</i>	✓	✓	✓
<i>FGFR4</i>	✓	✓	✓
<i>FH</i>	✓	-	-
<i>FLCN</i>	✓	-	-
<i>FLI1</i>	✓	-	✓
<i>FLT1</i>	✓	-	✓
<i>FLT3</i>	✓	-	✓
<i>FLT4</i>	✓	-	-
<i>FOXA1</i>	✓	-	-
<i>FOXL2</i>	✓	-	-
<i>FOXO1</i>	✓	-	-
<i>FOXP1</i>	✓	-	-
<i>FRS2</i>	✓	-	-
<i>FUBP1</i>	✓	-	-
<i>FYN</i>	✓	-	-
<i>GABRA6</i>	✓	-	-
<i>GATA1</i>	✓	-	-
<i>GATA2</i>	✓	-	-
<i>GATA3</i>	✓	-	-
<i>GATA4</i>	✓	-	-
<i>GATA6</i>	✓	-	-
<i>GEN1</i>	✓	-	-
<i>GID4</i>	✓	-	-
<i>GLI1</i>	✓	-	-
<i>GNA11</i>	✓	-	-
<i>GNA13</i>	✓	-	-
<i>GNAQ</i>	✓	-	-
<i>GNAS</i>	✓	-	-
<i>GPR124</i>	✓	-	-
<i>GPS2</i>	✓	-	-
<i>GREM1</i>	✓	-	-
<i>GRIN2A</i>	✓	-	-
<i>GRM3</i>	✓	-	-
<i>GSK3B</i>	✓	-	-
<i>H3F3A</i>	✓	-	-
<i>H3F3B</i>	✓	-	-
<i>H3F3C</i>	✓	-	-
<i>HGF</i>	✓	-	-
<i>HIST1H1C</i>	✓	-	-
<i>HIST1H2BD</i>	✓	-	-
<i>HIST1H3A</i>	✓	-	-
<i>HIST1H3B</i>	✓	-	-
<i>HIST1H3C</i>	✓	-	-

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
<i>HIST1H3D</i>	✓	-	-
<i>HIST1H3E</i>	✓	-	-
<i>HIST1H3F</i>	✓	-	-
<i>HIST1H3G</i>	✓	-	-
<i>HIST1H3H</i>	✓	-	-
<i>HIST1H3I</i>	✓	-	-
<i>HIST1H3J</i>	✓	-	-
<i>HIST2H3A</i>	✓	-	-
<i>HIST2H3C</i>	✓	-	-
<i>HIST2H3D</i>	✓	-	-
<i>HIST3H3</i>	✓	-	-
<i>HLA-A</i>	*	-	-
<i>HLA-B</i>	*	-	-
<i>HLA-C</i>	*	-	-
<i>HNF1A</i>	✓	-	-
<i>HNRNPK</i>	✓	-	-
<i>HOXB13</i>	✓	-	-
<i>HRAS</i>	✓	-	-
<i>HSD3B1</i>	✓	-	-
<i>HSP90AA1</i>	✓	-	-
<i>ICOSLG</i>	✓	-	-
<i>ID3</i>	✓	-	-
<i>IDH1</i>	✓	-	-
<i>IDH2</i>	✓	-	-
<i>IFNGR1</i>	✓	-	-
<i>IGF1</i>	✓	-	-
<i>IGF1R</i>	✓	-	-
<i>IGF2</i>	✓	-	-
<i>IKBKE</i>	✓	-	-
<i>IKZF1</i>	✓	-	-
<i>IL10</i>	✓	-	-
<i>IL7R</i>	✓	-	-
<i>INHA</i>	✓	-	-
<i>INHBA</i>	✓	-	-
<i>INPP4A</i>	✓	-	-
<i>INPP4B</i>	✓	-	-
<i>INSR</i>	✓	-	-
<i>IRF2</i>	✓	-	-
<i>IRF4</i>	✓	-	-
<i>IRS1</i>	✓	-	-
<i>IRS2</i>	✓	-	-
<i>JAK1</i>	✓	-	-
<i>JAK2</i>	✓	✓	✓

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
JAK3	✓	-	-
JUN	✓	-	-
KAT6A	✓	-	-
KDM5A	✓	-	-
KDM5C	✓	-	-
KDM6A	✓	-	-
KDR	✓	-	✓
KEAP1	✓	-	-
KEL	✓	-	-
KIF5B	✓	-	✓
KIT	✓	✓	✓
KLF4	✓	-	-
KLHL6	✓	-	-
KMT2B	*	-	-
KMT2C	*	-	-
KMT2D	*	-	-
KRAS	✓	✓	-
LAMP1	✓	✓	-
LATS1	✓	-	-
LATS2	✓	-	-
LMO1	✓	-	-
LRP1B	✓	-	-
LYN	✓	-	-
LZTR1	✓	-	-
MAGI2	✓	-	-
MALT1	✓	-	-
MAP2K1	✓	-	-
MAP2K2	✓	-	-
MAP2K4	✓	-	-
MAP3K1	✓	-	-
MAP3K13	✓	-	-
MAP3K14	✓	-	-
MAP3K4	✓	-	-
MAPK1	✓	-	-
MAPK3	✓	-	-
MAX	✓	-	-
MCL1	✓	-	-
MDC1	✓	-	-
MDM2	✓	✓	-
MDM4	✓	✓	-
MED12	✓	-	-
MEF2B	✓	-	-
MEN1	✓	-	-

## TruSight Oncology 500

---

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
<i>MET</i>	✓	✓	✓
<i>MGA</i>	✓	-	-
<i>MITF</i>	✓	-	-
<i>MLH1</i>	✓	-	-
<i>MLL</i>	✓	-	✓
<i>MLLT3</i>	✓	-	✓
<i>MPL</i>	✓	-	-
<i>MRE11A</i>	✓	-	-
<i>MSH2</i>	✓	-	✓
<i>MSH3</i>	✓	-	-
<i>MSH6</i>	✓	-	-
<i>MST1</i>	✓	-	-
<i>MST1R</i>	✓	-	-
<i>MTOR</i>	✓	-	-
<i>MUTYH</i>	✓	-	-
<i>MYB</i>	✓	-	-
<i>MYC</i>	✓	✓	✓
<i>MYCL1</i>	✓	✓	-
<i>MYCN</i>	✓	✓	-
<i>MYD88</i>	✓	-	-
<i>MYOD1</i>	✓	-	-
<i>NAB2</i>	✓	-	-
<i>NBN</i>	✓	-	-
<i>NCOA3</i>	✓	-	-
<i>NCOR1</i>	✓	-	-
<i>NEGR1</i>	✓	-	-
<i>NF1</i>	✓	-	-
<i>NF2</i>	✓	-	-
<i>NFE2L2</i>	✓	-	-
<i>NFKBIA</i>	✓	-	-
<i>NKX2-1</i>	✓	-	-
<i>NKX3-1</i>	✓	-	-
<i>NOTCH1</i>	✓	-	✓
<i>NOTCH2</i>	✓	-	✓
<i>NOTCH3</i>	✓	-	✓
<i>NOTCH4</i>	✓	-	-
<i>NPM1</i>	✓	-	-
<i>NRAS</i>	✓	✓	-
<i>NRG1</i>	✓	✓	✓
<i>NSD1</i>	✓	-	-
<i>NTRK1</i>	✓	-	✓
<i>NTRK2</i>	✓	-	✓
<i>NTRK3</i>	✓	-	✓



## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
<i>NUP93</i>	✓	-	-
<i>NUTM1</i>	✓	-	-
<i>PAK1</i>	✓	-	-
<i>PAK3</i>	✓	-	-
<i>PAK7</i>	✓	-	-
<i>PALB2</i>	✓	-	-
<i>PARK2</i>	✓	-	-
<i>PARP1</i>	✓	-	-
<i>PAX3</i>	✓	-	✓
<i>PAX5</i>	✓	-	-
<i>PAX7</i>	✓	-	✓
<i>PAX8</i>	✓	-	-
<i>PBRM1</i>	✓	-	-
<i>PDCD1</i>	✓	-	-
<i>PDCD1LG2</i>	✓	-	-
<i>PDGFRA</i>	✓	✓	✓
<i>PDGFRB</i>	✓	✓	✓
<i>PDK1</i>	✓	-	-
<i>PDPK1</i>	✓	-	-
<i>PGR</i>	✓	-	-
<i>PHF6</i>	✓	-	-
<i>PHOX2B</i>	✓	-	-
<i>PIK3C2B</i>	✓	-	-
<i>PIK3C2G</i>	✓	-	-
<i>PIK3C3</i>	✓	-	-
<i>PIK3CA</i>	✓	✓	✓
<i>PIK3CB</i>	✓	✓	-
<i>PIK3CD</i>	✓	-	-
<i>PIK3CG</i>	✓	-	-
<i>PIK3R1</i>	✓	-	-
<i>PIK3R2</i>	✓	-	-
<i>PIK3R3</i>	✓	-	-
<i>PIM1</i>	✓	-	-
<i>PLCG2</i>	✓	-	-
<i>PLK2</i>	✓	-	-
<i>PMAIP1</i>	✓	-	-
<i>PMS1</i>	✓	-	-
<i>PMS2</i>	✓	-	-
<i>PNRC1</i>	✓	-	-
<i>POLD1</i>	✓	-	-
<i>POLE</i>	✓	-	-
<i>PPARG</i>	✓	-	✓
<i>PPM1D</i>	✓	-	-

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
<i>PPP2R1A</i>	✓	-	-
<i>PPP2R2A</i>	✓	-	-
<i>PPP6C</i>	✓	-	-
<i>PRDM1</i>	✓	-	-
<i>PREX2</i>	✓	-	-
<i>PRKAR1A</i>	✓	-	-
<i>PRKCI</i>	✓	-	-
<i>PRKDC</i>	✓	-	-
<i>PRSS8</i>	✓	-	-
<i>PTCH1</i>	✓	-	-
<i>PTEN</i>	✓	✓	-
<i>PTPN11</i>	✓	-	-
<i>PTPRD</i>	✓	-	-
<i>PTPRS</i>	✓	-	-
<i>PTPRT</i>	✓	-	-
<i>QKI</i>	✓	-	-
<i>RAB35</i>	✓	-	-
<i>RAC1</i>	✓	-	-
<i>RAD21</i>	✓	-	-
<i>RAD50</i>	✓	-	-
<i>RAD51</i>	✓	-	-
<i>RAD51B</i>	✓	-	-
<i>RAD51C</i>	✓	-	-
<i>RAD51D</i>	✓	-	-
<i>RAD52</i>	✓	-	-
<i>RAD54L</i>	✓	-	-
<i>RAF1</i>	✓	✓	✓
<i>RANBP2</i>	✓	-	-
<i>RARA</i>	✓	-	-
<i>RASA1</i>	✓	-	-
<i>RB1</i>	✓	-	-
<i>RBM10</i>	✓	-	-
<i>RECQL4</i>	✓	-	-
<i>REL</i>	✓	-	-
<i>RET</i>	✓	✓	✓
<i>RFWD2</i>	✓	-	-
<i>RHEB</i>	✓	-	-
<i>RHOA</i>	✓	-	-
<i>RICTOR</i>	✓	✓	-
<i>RIT1</i>	✓	-	-
<i>RNF43</i>	✓	-	-
<i>ROS1</i>	✓	-	✓
<i>RPS6KA4</i>	✓	-	-

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
<i>RPS6KB1</i>	✓	✓	✓
<i>RPS6KB2</i>	✓	-	-
<i>RPTOR</i>	✓	-	-
<i>RUNX1</i>	✓	-	-
<i>RUNX1T1</i>	✓	-	-
<i>RYBP</i>	✓	-	-
<i>SDHA</i>	✓	-	-
<i>SDHAF2</i>	✓	-	-
<i>SDHB</i>	✓	-	-
<i>SDHC</i>	✓	-	-
<i>SDHD</i>	✓	-	-
<i>SETBP1</i>	✓	-	-
<i>SETD2</i>	✓	-	-
<i>SF3B1</i>	✓	-	-
<i>SH2B3</i>	✓	-	-
<i>SH2D1A</i>	✓	-	-
<i>SHQ1</i>	✓	-	-
<i>SLIT2</i>	✓	-	-
<i>SLX4</i>	✓	-	-
<i>SMAD2</i>	✓	-	-
<i>SMAD3</i>	✓	-	-
<i>SMAD4</i>	✓	-	-
<i>SMARCA4</i>	✓	-	-
<i>SMARCB1</i>	✓	-	-
<i>SMARCD1</i>	✓	-	-
<i>SMC1A</i>	✓	-	-
<i>SMC3</i>	✓	-	-
<i>SMO</i>	✓	-	-
<i>SNCAIP</i>	✓	-	-
<i>SOCS1</i>	✓	-	-
<i>SOX10</i>	✓	-	-
<i>SOX17</i>	✓	-	-
<i>SOX2</i>	✓	-	-
<i>SOX9</i>	✓	-	-
<i>SPEN</i>	✓	-	-
<i>SPOP</i>	✓	-	-
<i>SPTA1</i>	✓	-	-
<i>SRC</i>	✓	-	-
<i>SRSF2</i>	✓	-	-
<i>STAG1</i>	✓	-	-
<i>STAG2</i>	✓	-	-
<i>STAT3</i>	✓	-	-
<i>STAT4</i>	✓	-	-

## TruSight Oncology 500

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
STAT5A	✓	-	-
STAT5B	✓	-	-
STK11	✓	-	-
STK40	✓	-	-
SUFU	✓	-	-
SUZ12	✓	-	-
SYK	✓	-	-
TAF1	✓	-	-
TBX3	✓	-	-
TCEB1	✓	-	-
TCF3	✓	-	-
TCF7L2	✓	-	-
TERC	✓	-	-
TERT	✓	-	-
TET1	✓	-	-
TET2	✓	-	-
TFE3	✓	-	-
TFRC	✓	✓	-
TGFBR1	✓	-	-
TGFBR2	✓	-	-
TMEM127	✓	-	-
TMPRSS2	✓	-	✓
TNFAIP3	✓	-	-
TNFRSF14	✓	-	-
TOP1	✓	-	-
TOP2A	✓	-	-
TP53	✓	-	-
TP63	✓	-	-
TRAF2	✓	-	-
TRAF7	✓	-	-
TSC1	✓	-	-
TSC2	✓	-	-
TSHR	✓	-	-
U2AF1	✓	-	-
VEGFA	✓	-	-
VHL	✓	-	-
VTCN1	✓	-	-
WISP3	✓	-	-
WT1	✓	-	-
XIAP	✓	-	-
XPO1	✓	-	-
XRCC2	✓	-	-
YAP1	✓	-	-

## TruSight Oncology 500

---

Gene symbol	Small variants	Focal amplifications	Fusions (from RNA)
YES1	✓	-	-
ZBTB2	✓	-	-
ZBTB7A	✓	-	-
ZFHX3	✓	-	-
ZNF217	✓	-	-
ZNF703	✓	-	-
ZRSR2	✓	-	-

\* Small variants found in gVCF file only

*Some of the genomic variants, genes, nucleic acid sequences, or genomic regions on this list, and their use in specific applications, may be protected by patents. Customers are advised to determine whether they are required to obtain licenses from the party that owns or controls such patents in order to use the product in Customer's specific application. Unless expressly stated otherwise in writing by Illumina, Customer, and not Illumina, is solely responsible for obtaining such licenses.*